




**Chase  
Research  
Cryogenics**

**World leaders in  
sub-Kelvin cryogenics**

## **CONTINUOUS $^4\text{He}$ SORPTION COOLER**

### **TYPE CC4**

## **GENERIC INSTALLATION AND OPERATING INSTRUCTIONS**



**Photo shows a typical CRC CC4 sorption cooler**

**CRC Ltd, Cool Works, Unit 2 Neepsend Ind. Est., Parkwood Rd, Sheffield S3 8AG, UK.**

**Tel: +44 114 278 0711**

**Director: Dr. S.T. Chase**

**Registered in England & Wales, No. 4643351**

**Secretary: Dr. L.C. Kenny**

**VAT registration No. GB 763 8558 84**

## Contents

1. GENERAL HANDLING.....	3
2. SAFETY OF CHASE RESEARCH CRYOGENICS PRODUCTS.....	3
2.1. Pressure Equipment Directive 97/23/EC (Pressure Equipment Regulations 1999).....	3
2.2. Pressure Systems Safety Regulations 2000.....	3
2.3. Safe Operation.....	4
2.4. Risk Assessment .....	4
3. A BRIEF DESCRIPTION OF THE SORPTION COOLER UNIT.....	6
4. INSTALLATION.....	7
4.1. Mechanical .....	7
4.3. Electrical .....	8
5. ATTACHING YOUR EXPERIMENT TO THE SORPTION COOLER.....	9
5.1. Radiation shielding .....	9
6. OPERATION: quick-start guide.....	10
6.1. Cooldown and commence continuous operation.....	10
7. OPTIMISING THE PERFORMANCE OF YOUR SORPTION COOLER .....	11
7.1. Varying the cycle time .....	11
7.2. Varying the cycle steps.....	12
7.3. PID stabilisation of the CC4 .....	13
9. STANDARD Pin-Out Assignments .....	14
10. cc4 OPERATIONAL SEQUENCE FOR INITIAL COOLDOWN AND EXTENDED OPERATION .....	15

**THIS GENERIC OPERATING MANUAL** describes how to install and operate a CRC CC4 sorption cooler. It is accompanied by an Excel file that contains the validation test data and the calibration files that are **specific** to the sorption cooler unit that you have purchased.

You are advised to make a note of the location of the Excel file specific to your sorption cooler unit. CRC can provide another copy on request if you send us the Serial Number engraved around the rim of your cooler's main plate.

This manual was last updated in January 2022.

## 1. GENERAL HANDLING

### **WARNING!**

#### **CRC SORPTION COOLERS CONTAIN HELIUM GAS AT HIGH PRESSURE.**

**Do not crush, twist or bend the unit. Avoid applying mechanical stresses. Do not heat the unit above room temperature. Keep in a sealed cryostat, or in the shipping box and brace in which it came.**

**Do not hold or lift the unit by means of the cold heads.**

**Do not tamper with the copper capillary fill tubes.**

**Avoid the use of acid fluxes when soldering in the vicinity of the cooler. Chloride based fluxes will corrode stainless steel and could damage your cooler.**

**Always allow the unit to warm to room temperature before allowing air into your cryostat. Opening the air valve of a cold cryostat is likely to damage the unit.**

After unpacking the sorption cooler according to the instructions supplied, the sorption cooler should be immediately transferred into the host cryostat. The shipping brace doubles as a stand for the sorption cooler, though when used as a stand, the screws through the top aluminium plate into the cold heads should NOT be in place. When picking the sorption cooler up, it should be firmly grasped by the main plate.

## 2. SAFETY OF CHASE RESEARCH CRYOGENICS PRODUCTS

### **2.1. Pressure Equipment Directive 97/23/EC (Pressure Equipment Regulations 1999)**

This CRC sorption cooler unit is manufactured in accordance with Sound Engineering Practice. The volume and gas pressure within the sorption cooler are such that the equipment falls below the lower classification limit in Annex II of the Pressure Equipment Directive. Hence the requirements for Conformity Assessment do not apply and no Declaration of Conformity can be made, or CE marking applied.

The sorption cooler is covered by Article 3 Paragraph 3 of the Pressure Equipment Directive, which states: "Pressure equipment and/or assemblies below or equal to the limits in sections 1.1, 1.2 and 1.3 and section 2 respectively must be designed and manufactured in accordance with the sound engineering practice of a Member State in order to ensure safe use. Pressure equipment and/or assemblies must be accompanied by adequate instructions for use and must bear markings to permit identification of the manufacturer or of his authorized representative established within the Community. Such equipment and/or assemblies must not bear the CE marking referred to in Article 15."

### **2.2. Pressure Systems Safety Regulations 2000**

This sorption cooler unit does not contain a pressure x volume product exceeding 250 bar-litres hence PSSR regulations 5(4), 8-10 and 14 do not apply. This means that the system does not require a written scheme of examination. The sorption cooler is not 'mobile' in the sense intended in the PSSR hence *the owner* has duties under these regulations to ensure that a) the safe operating limits are not exceeded; b) the unit is operated in accordance with these instructions; c) the unit is

returned to Chase Research Cryogenics Ltd in the event that any maintenance is required. The sorption cooler contains no user-serviceable parts.

### **2.3. Safe Operation**

The safe operating temperature range of this sorption cooler is 0 to 320 K.

### **2.4. Risk Assessment**

CRC sorption coolers contain Helium gas under pressure. The stored energy of the system is less than 50 bar litres. All system components are integrity tested during manufacture; the slightest leak will make the sorption cooler lose its stored gas and cease to function. A unit that has leaked presents no risks whatever to the user; the following risk assessment applies therefore only to functional units.

#### ***Hazards and consequences***

Accidental damage to the sorption cooler unit could result in the sudden release of pressurised gases, causing mechanical failure of the unit and potential injury (or damage to surrounding instruments) from ejected debris.

Possible events leading to failure are: overheating of the unit, for example in a fire; dropping or crushing of the unit; twisting or bending of the gas tubes. Mechanical damage to the unit is most likely to occur during assembly of the instrument of which the sorption cooler forms part.

#### ***Risks without controls in place***

It is extremely unlikely that the above events will lead to danger. Chase Research Cryogenics Ltd has produced more than two hundred sorption cooler units of various designs, which are in use for a range of applications worldwide. To date there has never been a sudden failure of a sorption cooler unit – indicating that with normal use (including inevitable handling mishaps) the units have an excellent safety record. User experience to date shows that accidental mechanical damage to sorption cooler units is likely to result in slow leaks, not sudden failures.

#### ***Controls in place***

The controls that are in place to eliminate (as far as reasonably practicable) the risks arising from mechanical damage to a sorption cooler unit are:

- This written instruction manual, containing warnings about the potential risks arising from damage to the unit and alerting the user to more risky operations;
- Instructions that the unit should not be used if it has been subjected to overheating, dropping, crushing, bending or twisting;
- A warning label on the transit box that the instructions should be read prior to handling the unit.

The applications for which sorption cooler units are intended make it impossible to place warning labels on the unit itself. However if the sorption cooler is incorporated into another instrument, that instrument should carry a warning label to alert the user that the sorption cooler contains no user-serviceable parts and should not be disassembled.

#### ***Risks with controls in place***

Providing users read and follow this instruction manual the risks are negligible.

## 2.5. Safe Transportation

The unit must be correctly re-installed in its shipping brace before transportation either by road or air. Follow the unpacking instructions provided in reverse order and contact CRC for advice if needed. The cooler+brace should be securely packed into the cooler's rigid shipping box, which contains polystyrene or similar cushioning material to firmly hold the unit and prevent any movement within the box. Add extra pieces of packing foam to ensure that the heavy pumps cannot move in transit. Place the rigid shipping box into a much larger outer carton surrounded by at least 5-10 cm soft fill on all sides, this should ensure that if dropped the impact shock will be absorbed by the overpack, and not transmitted to the cooler itself.


For transportation purposes CRC sorption coolers are classified under UN 2857, Class 6A (6=Other articles containing gas under pressure; A=asphyxiant gas). Special Provision 119 applies: *Refrigerating machines include machines or other appliances which have been designed for the specific purpose of keeping food or other items at a low temperature in an internal compartment, and air conditioning units. Refrigerating machines and refrigerating machine components are not subject to the provisions of ADR if they contain less than 12 kg of gas in Class 2, group A or O according to [2.2.2.1.3](#), or if they contain less than 12 litres ammonia solution (UN No. 2672).*

This special provision means that the cooler does not need to be either labelled or shipped as hazardous or dangerous goods, whether transported by road or by air, provided that it is appropriately packed.

Ensure that the outer carton is clearly labelled as fragile with 'This Way Up' labels and add Shock Indicators if you are handing the consignment to third-party couriers, as these should deter rough handling.

### 3. A BRIEF DESCRIPTION OF THE SORPTION COOLER UNIT

The CC4 sorption cooler is shown in the figure below. Continuous cooling of the 1K cold head is achieved by alternately cycling the two  $^4\text{He}$  module at intervals, such that the 1K cold head is always being cooled by one of the modules. In use, the sorption cooler is inverted, i.e. the cold head will be at the bottom. The main plate needs to be thermally sunk to the cold head of a cryocooler at 4K or below, see section 4.1 for more information. The 1K-head is used to extract heat from the experiment, see section 5 for more information. The pumps and heat switches can reach up to 50K during operation, these need to be radiation shielded from the head, for more information see section 5.1. All electrical connections are mounted on the main plate on a stack of MDM connectors, see section 4.2 for more information.



The following short names for the various parts of the sorption cooler are used throughout this user manual:

Short name used in this manual	Refers to the sorption cooler part
1K-head	The cold head of the split condenser module
4-head	$^4\text{He}$ module head
4-Pump	$^4\text{He}$ module pump
4-Switch	Heat Switch for the $^4\text{He}$ pump

The two  $^4\text{He}$  modules (which are identical) will be referred to as Module A and Module B.

## 4. INSTALLATION

### 4.1. Mechanical

**Before installing the unit in your cryostat, be sure to carefully remove any pieces of foam packing from around the pumps, as mentioned in the unpacking instructions.**

**There should be no need to touch the heat switches or heat straps during installation or normal operation of the sorption cooler. The heat switches can be easily damaged, and if bent or twisted are likely to fail.**

This sorption cooler is designed to be pre-cooled using a mechanical cryocooler such as a GM or pulse tube operating at 4K or below. The thermal link to the cryocooler should be made from gold plated copper to ensure excellent thermal contact between the sorption cooler and the cryocooler. Mounting holes are provided on the main plate for attaching the sorption cooler to your cryostat cold plate. To attach the sorption cooler to the 4K stage of the cryocooler there are twelve M4 clearance holes distributed around the periphery of the circular main plate (note: UNC6 clearance holes can be substituted at customer request). A CAD file of your sorption cooler can be provided on request to aid with interfacing.

**Because the cooling down of the heads depends upon gas convection, and on liquid helium collecting in the heads fed by gravity, the sorption cooler *must* be kept close to vertical with the heads downwards.**



**Ensure spring washers are under every bolt head, these will take out differential thermal contraction that might otherwise cause loosening of the bolts, and thus compromise thermal contact.**

### 4.3. Electrical

All electrical connections are brought out to stack of three MDM-SSP connectors (2 x 15 pin, 1 x 9 pin) mounted onto the main plate. Pin-outs are listed at the end of this instruction manual. The table below summarises the temperature sensors installed on the CC4.

ITEM	Calibration	Options
<b>1K-head RuO<sub>2</sub></b>	Generic Lakeshore Cryotronics or Scientific Instruments Inc.	Individually calibrated sensors available on request
<b>4-head RuO<sub>2</sub></b>	Generic housekeeping sensors, calibration supplied by CRC	Lakeshore sensors available on request
<b>Switch diode</b>	Generic – supplied by CRC Ltd	
<b>Pump diode</b>	Generic – supplied by CRC Ltd	

Voltage / current requirements for driving the heater and thermometers are summarised in the table below:

ITEM	NUMBER	IMPEDANCE/ JUNCTION VOLTAGE	VOLTAGE/ CURRENT
<b>4-pump heater</b>	1	300 $\Omega$ approx.	20 to 25 V
<b>Heat switch heater</b>	1	10k $\Omega$	5 V
<b>Diode thermometers</b>	3 or 4	0.5 to 1.8V	10 $\mu$ A DC
<b>4-head RuO<sub>2</sub> thermometer</b>	1	1k $\Omega$ to 3k $\Omega$	1 $\mu$ A max.

Generic (i.e. standard calibration) RuO<sub>2</sub> sensors from Lakeshore Cryotronics or Scientific Instruments Inc. are the default option on the head of all CRC sorption coolers. Individually calibrated ‘CERNOX’ or RuO<sub>2</sub> sensors are only fitted (at additional cost) at the customer’s request. The thermometer on the 1K-head is operated as a 4-wire device and should ideally be driven by an AC current no greater than 1 $\mu$ A. Calibration tables for all thermometer sensors are in the Excel data file that accompanies each unit. Generic diode calibration curves for the pump diode and heat switch diode are supplied as standard. The diode thermometers require excitation with currents of 10 $\mu$ A DC.


The heat switch heater typically requires about 5V to keep the switch in the ‘ON’ state with the absorber pod at greater than  $\sim$ 20K, and it will cool to the off state ( $T < 10$  K) in ten to fifteen minutes. The pump heater impedance is typically 300 $\Omega$ . To operate the sorption cooler it is necessary to periodically warm the 4-pumps to between 40 and 45k. A heater current of up to 100 to 130mA will heat the 4-pump rapidly; lower heater currents will result in slower heating. Try to ensure that the lead-in wiring to the pump heater is not unduly dissipative and is appropriately heat-sunk.

## 5. ATTACHING YOUR EXPERIMENT TO THE SORPTION COOLER.

This model of sorption cooler provides just one point at which heat may be extracted from a user's experiment mounted on a separate cold table, which is the 1K-head. There are several M3 holes tapped in this surface for thermal connections between your experiment and the sorption cooler (note UNC4 holes can be substituted at customer request).

**While fixing experimental equipment to the cold head, extreme care should be taken not to torque or bend the gas pipes. Always support the cold head against the applied torque.**

Under no load, and with the main plate at  $\sim 4\text{K}$ , the cold head will run at an average temperature of about  $900\text{mK}$ . The temperature fluctuates periodically due to the antiphase cycling of the 4-modules; the average temperature and the amplitude of the fluctuations will depend on the software used to control the CC4 and there are different optimisation strategies that are possible, depending on whether low temperature or temperature stability is more important. This is discussed further in Section 7. When the 1K-head is loaded the temperature increases, typical data are given in Figure 1 below. Load data for your specific sorption cooler are supplied in the Excel data file that accompanies your unit.



### 5.1. Radiation shielding

The 1K-head, and any cold table/experimental equipment/detector assembly you attach, must be properly radiation shielded at around  $4\text{K}$  to achieve satisfactory operation. The cold table should be thermally connected to the 1K-head with a copper heat strap. No other mechanical attachments to the sorption cooler unit are necessary for satisfactory operation. If your sorption cooler performance is not meeting the specification, this is likely to be due to a radiation load. Check your radiation shielding and consider adding extra multi-layered insulation around your radiation shields, or around the pump side ("warm side") of the CC4.

## 6. OPERATION: QUICK-START GUIDE

### 6.1. Cooldown and commence continuous operation

The basic cooldown steps for the CC4 are as follows:

- Turn on the GM cryocooler and wait until both heat switches have turned off.
- Turn on the pump heaters to keep both pumps warm (around 45K) until both 4-heads and the 1K-head are cold (4K or below)
- Turn the heat off pump A, turn on heat switch A to cool 4-head A to ~1K. The 1K-head will also cool to ~1K at this point.
- 15 minutes later: Turn the heat off pump B, turn on heat switch B to cool 4-head B to below 1K.
- With all heads now below 1K the periodic re-cycling can begin.

This cooldown sequence is illustrated in Figure 2 below, which shows a typical cooldown from room temperature to 1K (at ~1350 minutes). The key steps to observe are:

- Heat switches turn off just after  $t=1200$
- Pump heaters are turned on shortly after this
- All heads are at ~4K around  $t=1340$
- Switch A turned on just before  $t=1350$  for final cooldown of 4-head A and 1K-head
- Switch B turned on 15 minutes after for final cooldown of 4-head B
- Ready to begin cycling at  $t=1350$




Figure 2

## 6.2 Established alternate cycling for continuous operation.

The cycle time that is optimum for your CC4 will depend on your cryocooler's power and the heat load that your experiment imposes, but is likely to be in the range 30-100 minutes. The examples in this manual use 100 minutes cycle time.

The cycling steps are explained more fully in Section 8 and can be summarised as follows:

- Switch A off
- Heat pump A (the 1K-head is being cooled by 4-head B)
- Keep pump A warm while 4-head A cools
- Switch A on
- Switch B off
- Heat pump B (the 1K-head is being cooled by 4-head A)
- Keep pump B warm while 4-head B cools
- Switch B on

This sequence can be repeated for as long as required.

## 7. OPTIMISING THE PERFORMANCE OF YOUR SORPTION COOLER

Running the CC4 with programmable power supplies under software control enables a range of options to tailor the CC4 performance to the user's own set-up and application. Options include:

- Varying the cycle time;
- Varying the timing of individual cycle steps;
- Varying the rates at which pumps and/or heat switches are heated.

You are recommended to experiment with your own system to optimise its performance in your setup, for your application.


### 7.1. Varying the cycle time

The cycle time cannot be reduced below the time taken to heat and cool the pumps and the heat switches. Very rapid heating cycles are only possible if the CC4 is run using a powerful GM/PT cryocooler that is able to rapidly remove the heat dumped by the pump when the heat switch is turned on. This heat pulse temporarily increases the temperature of the main plate, and if it gets too warm (above ~10K) there is a danger that the other heat switch (the 'off' switch) will also turn on. Disrupting the cycle. In general, for better temperature stability and low running temperature, we advise gentle turn-on of pump and heat switch heaters (see below).

The cycle time cannot be lengthened beyond the expiry time of either He4 module in the CC4, i.e. the time taken to evaporate all of the Helium 4 liquid in the 4-head. When a module expires it will start to warm up, and will warm the 1K-head before the other module takes over the cooling. The expiry time will depend on the load imposed on the head, both by parasitic loads internal to the CC4, and by external loads due to the user's experiment. Your CC4's technical datafile will specify a recommended cycle time for your specific unit, when working under a specified applied external load. If you find that your experiment causes the modules to expire, you should seek to reduce the load (i.e. improve radiation shielding and heat sinking), and then try reducing the cycle time.


## 7.2. Varying the cycle steps

Figure 3 shows an example of a cycling sequence (100 minute cycle time) optimised to give a low base temperature for the CC4. In this example the pump heaters are turned on gently at first, and they are turned off as soon as the 4-heads stop cooling rapidly. This minimises the heat loading on the CC4.



**Figure 3**

In the second example shown in Figure 4 the cycling sequence is optimised to reduce temperature fluctuations at the 1K-head.



**Figure 4**

The cycling time in Figure 4 is again 100 minutes, but here, the heat switches are turned on gently at first, and the pumps are kept warm for a longer portion of each cycle. This means that the heat loading on the CC4 is greater overall, however the rate of heat transfer to the main plate when the cycle is switched is reduced, meaning that the temperature fluctuations transmitted to the 1K-head are also reduced.

### **7.3. PID stabilisation of the CC4**

Finally, the CC4 wiring is configured to enable PID stabilisation of the temperature through the application of a varying external load via a feedback loop to compensate the inbuilt fluctuations. Many different PID optimisation schemes are possible and the user is encouraged to experiment to find the one that best suits their own set-up.

PID optimisation does require the application of loads to the 1K-head and so the CC4 will run at a higher, though more stable, temperature when it is used.

## 9. STANDARD PIN-OUT ASSIGNMENTS

### He4 modules

Function		female 15 pin	Drive current or voltage
	TYPE	PIN #	
RTD THERMOMETER V+	generic 6951	1	100nA AC
RTD THERMOMETER V-		9	Or low voltage
RTD THERMOMETER I+		2	Driver
RTD THERMOMETER I-		10	e.g. V<0.5mV
not used		3	
PUMP DIODE I+	DC2018	4	10 $\mu$ A
PUMP DIODE I-		11	
SWITCH DIODE I+	DC2018	5	10 $\mu$ A
SWITCH DIODE I-		12	
not used		6	
not used		13	
SWITCH HEATER I+	10kOhm	7	2.5 to 3 Volts
SWITCH HEATER I-		14	
PUMP HEATER I+	300 Ohm	8	50 to 100 mA
PUMP HEATER I-		15	13 to 15 V

Function		female 9 pin	Drive current or voltage
	TYPE	PIN #	
RTD THERMOMETER V+	Lakeshore or Scientific Instruments	1	100nA AC
RTD THERMOMETER V-		6	Or low voltage
RTD THERMOMETER I+		2	Driver
RTD THERMOMETER I-		7	e.g. V<0.5mV
nc		3	
nc		4	
nc		8	
PID heater	100 kOhm	5	50 to 100 mA
		9	0 to 5 V

Resistance thermometer

Diode thermometer

Low power heater (a few mW)

High power heater (up to about  
2W)



## 10. CC4 OPERATIONAL SEQUENCE FOR INITIAL COOLDOWN AND EXTENDED OPERATION

Step	Pump A heater	Pump B heater	A Heat switch	B Heat switch	Working status	Notes
1	OFF	OFF	OFF	OFF	Wait for both heat switches to turn off	Pre-cooling with GM/PT
2	ON	ON	OFF	OFF	Wait until both 4-head temperatures are less than 4K	Heat both pumps to ~45 K, turn power down (or off) to keep them at 40-45K
3	OFF	OFF	ON	OFF	Turn on heat switch A, 4-head A and 1K-head cool down	heads should cool to ~1K fairly rapidly
4	OFF	OFF	ON	ON	Turn on heat switch B, 4-head B cools down	All heads should now be at ~1K
<b>Begin alternate cycling for extended operation. Repeat steps 5 to 11 at your chosen cycle interval</b>						
5	OFF	OFF	OFF	ON	Turn off heat switch A, 4-head A temperature rises to ~4K	The 1K-head is being cooled by Module B
6	ON	OFF	OFF	ON	Turn on Pump A and Heat to 45K, 4 head A desorption	A more gradual heat will reduce temperature fluctuations but raise the average temperature
7	OFF	OFF	OFF	ON	Turn off Pump A heater, Wait for thermal relaxation time	The time taken for this step will depend on the cryocooler power. Typically ~10 minutes
8	OFF	OFF	ON	ON	Turn on heat switch A, 4 head A starts to cool down	A gradual heat switch turn on will reduce temperature fluctuations
9	OFF	OFF	ON	OFF	Turn off heat switch B, 4head B temperature rises to~ 4K	The 1K-head is being cooled by Module A
10	OFF	ON	ON	OFF	Turn on Pump B and Heat to 45K, 4head B desorption	See step 6
11	OFF	OFF	ON	OFF	Turn off Pump B heater, Wait for thermal relaxation time	See step 7